


“Minicircle DNA will become a cornerstone of modern gene and cell therapy“

31 October 2025 | Views | By Ankit Kankar

PlasmidFactory GmbH, based in Bielefeld, Germany, is recognised as a leading CDMO for the production of plasmid and Minicircle DNA. Over the past 25 years, the company has evolved from a small academic spin-off into a global reference for DNA manufacturing in gene and cell therapy. Following the launch of its new GMP facility and the clinical success of several Minicircle-based therapies, PlasmidFactory is entering a new phase—expanding its impact from plasmid DNA manufacturing toward next-generation, virus-free gene delivery. In an interaction with BioSpectrum India Dr Ram Shankar, Chief Scientific Officer at PlasmidFactory, shares his insights into the company’s latest developments and the growing importance of Minicircle technology in cell and gene therapy.



You recently presented data on Minicircle-based CAR-T cells at the CAR-TCR Summit in Boston. What makes Minicircles such a breakthrough for gene and cell therapy?

Minicircles represent the logical next step in non-viral gene transfer. Traditional plasmids contain bacterial backbone sequences—origins of replication, antibiotic resistance genes, and other elements that are not only unnecessary for therapy but can also compromise safety and performance. Minicircles remove these bacterial sequences entirely, leaving only the therapeutic expression cassette in a small, supercoiled DNA molecule.

This design has profound advantages: higher transfection efficiency, stronger and more stable expression, and significantly lower cytotoxicity. Because there are no bacterial genes, the risk of immune activation or gene silencing is minimised. The result is a cleaner, safer, and more efficient DNA vector—ideally suited for advanced therapies such as CAR-T, TCR-T, or TIL approaches.

How do Minicircles address some of the limitations of viral vectors, such as lentivirus?

Viral vectors have enabled the first wave of gene therapies, but they come with several drawbacks—complex production, high cost, limited payload capacity, and concerns around insertional mutagenesis. Minicircles, in combination with transposon systems like *Sleeping Beauty* or *PiggyBac*, offer a virus-free alternative that is scalable, cost-effective, and safer.

In recent preclinical and clinical studies, we have shown that CAR-T cells engineered with Minicircle DNA and *Sleeping Beauty* transposase mRNA can eradicate tumors as efficiently as viral-vector CAR-T cells, while integrating preferentially into genomic safe harbors. This reduces the risk of insertional mutagenesis by more than sixfold compared to lentiviral vectors. Moreover, the production of Minicircles is much faster and more economical—there is no need for complex viral packaging or biosafety-level infrastructure.

Several clinical programmes already use Minicircles. Could you share some examples?

Yes, Minicircles produced at PlasmidFactory have already entered multiple clinical trials, including *CARAMBA*, *TranspoCART19*, and *LION-1*. These programmes explore virus-free CAR-T cell generation using transposon systems and have demonstrated strong efficacy and favorable safety profiles.

Our technology has therefore moved well beyond proof-of-concept. It is clinically validated and now ready for late-stage development and commercial translation. We are proud that Minicircle DNA produced by PlasmidFactory contributes directly to these pioneering studies, demonstrating that virus-free engineering is not just possible—it is the future.

What role does manufacturing quality play in enabling this transition from research to clinic?

Manufacturing quality is absolutely central to what we do. For 25 years, PlasmidFactory has set the benchmark for high-purity DNA. We offer three quality grades—Research, High Quality (HQ), and GMP—to ensure a seamless transition from early discovery to clinical and commercial production.

Our new GMP facility in Bielefeld, which will be fully operational in 2024, allows us to produce both plasmid and Minicircle DNA at large scale using single-use equipment and rigorous quality controls. Each product is made under strict segregation—only one plasmid or Minicircle is produced per suite—to avoid any risk of cross-contamination.

Our proprietary purification process delivers DNA with extremely high supercoiled (ccc) content and minimal impurities. These standards meet or exceed current regulatory expectations for advanced therapy starting materials.

Your Minicircle technology is patented. How does it differ from other non-viral systems available today?

The uniqueness lies both in the molecular design and in the process. Our patented system produces Minicircles by intramolecular recombination from a well-characterised parental plasmid. During the process, the bacterial backbone forms a miniplasmid that is efficiently removed during purification, leaving only the therapeutic Minicircle DNA.

The method is robust, scalable, and license-free for our clients, which is an important advantage. It also allows the production of complex constructs without cargo limitations—supporting polycistronic designs or regulatory elements for fine-tuned expression. Because the same platform is compatible with both research-grade and GMP manufacturing, our customers can easily scale their programmes without changing the vector backbone or process.

Beyond CAR-T, where else do you see Minicircle DNA having impact?

The potential applications are vast. Besides CAR-T and TCR-T cells, we see strong interest in AAV production and RNA manufacturing, where Minicircle DNA serves as a cleaner, more efficient template compared to conventional plasmids. It eliminates unwanted bacterial sequences, improving both safety and productivity of viral vector and mRNA production processes.

Beyond these, Minicircles are also gaining attention in natural killer (NK) cell engineering, induced pluripotent stem cells, and in vivo gene delivery. A particularly exciting frontier is direct in vivo immune-cell engineering, where Minicircles can be formulated with lipid or polymer nanoparticles (LNPs/PNPs) for systemic, non-viral gene transfer. This could one day enable in vivo generation of CAR-T or TIL cells, avoiding the need for complex ex vivo manufacturing entirely.

What trends do you observe in the broader CDMO landscape for DNA manufacturing?

We are witnessing an evolution from conventional plasmid production toward specialised DNA formats tailored for advanced therapies—Minicircles being a prime example. Sponsors now require CDMOs that not only manufacture DNA but also understand the underlying biology and regulatory context.

At PlasmidFactory, we combine 25 years of manufacturing experience with deep scientific expertise. Our scientists work closely with clients from vector design through analytical characterisation, ensuring that the DNA meets both functional and regulatory requirements. With over 3,500 constructs produced and a 99.9 per cent success rate, we offer a level of reliability and know-how that few others can match.

You have recently strengthened your global partnerships. How important is international collaboration, particularly in Asia and India?

It is very important. India has a rapidly growing biotechnology and vaccine industry with strong interest in RNA and DNA technologies. Together with our local partner, Dr Nagaraj Rao of RRR Labs, Navi Mumbai, we are building bridges between German DNA manufacturing expertise and India's dynamic biopharma ecosystem.

This partnership supports Indian researchers and developers by giving them access to the same high-quality DNA—both plasmid and Minicircle—that leading global companies use for their advanced therapy programs. We see great potential for collaboration in the areas of CAR-T development, RNA therapeutics, and DNA vaccine research.

From a scientific point of view, what are your priorities for R&D?

Our R&D covers the entire DNA production chain – from strain and construct optimisation to upstream bioprocessing, purification, and analytical control. We focus on process innovation, plasmid stability, and quality assurance, continuously refining fermentation, purification, and single-use technologies to maximise yield and purity. A key area is the development of next-generation vectors, including non-viral and resistance-gene-free systems such as Minicircles, and their application in AAV and RNA production.

Finally, where do you see the future of PlasmidFactory and of Minicircle technology in the next five years?

We believe that Minicircle DNA will become a cornerstone of modern gene and cell therapy. As the industry moves away from viral systems, demand for safe, efficient, and scalable non-viral vectors will grow exponentially. PlasmidFactory is ready for that future—with proven technology, GMP capacity, and a passionate team dedicated to advancing the science of DNA manufacturing.

In the next five years, we expect to see Minicircles not only in ex-vivo CAR-T programmes but also in *in vivo* applications using LNP-based delivery. This will truly transform how we think about gene therapy—simpler, safer, and accessible to more patients worldwide.

Ankit Kankar

ankit.kankar@mmactiv.com